

## pluriStrainer® 100 µm (filtros de células)



	Order No.	Sterility	Packaging	Variant
<input type="checkbox"/>	43-50100-51	Estéril	Individual single packed	50 pcs, sterile, PET-Mesh
<input type="checkbox"/>	43-50100-50	No estéril	50 pcs. / Bag	50 pcs, non-sterile, PET-Mesh
<input type="checkbox"/>	43-50100-03	Estéril	Individual single packed	25 pcs, sterile, PET-Mesh
<input type="checkbox"/>	43-57100-50	No estéril	50 pcs. / Bag	50 pcs, non-sterile, Nylon-Mesh
<input type="checkbox"/>	43-57100-03	Estéril	Individual single packed	25 pcs, sterile, Nylon-Mesh
<input type="checkbox"/>	43-57100-51	Estéril	Individual single packed	50 pcs, sterile, Nylon-Mesh

### Descripción

The pluriStrainer® is a sieving device (or cell strainer) for a wide range of laboratory applications where the filtration and the purification of liquids is required. Its unique design features improved ventilation to avoid clogging for a smart sample preparation. pluriStrainer® are stackable to allow direct filtration with different mesh sizes and also can be inverted to recover the sieved material. For large sample volumes it can be combined with a funnel. It fits all major brand 50 mL centrifuge (conical) tubes. The combination of the pluriStrainer® and the Connector Ring allow using low pressure which supports the filtration. Furthermore it is possible to block the flow through the sieve and hold back the liquid on top of the mesh. The liquid on top of the mesh can be used for sample preparation and incubation such as tissue dissociation or the treatment of cells with cell lysis reagents (e.g. Buffer RLT, TRIzol®).

### Aplicación

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Enrichment of specific cells using pluriBead technology
- Preparation of a single cell suspension of blood cells from bone marrow, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry (FACS™)
- Faster and easier alternative to gauze filtration
- Sieving of complex and viscous liquids (in combination with Connector Ring & Syringe for vacuum)
- Short time incubation of cells / cytokine (in combination with pluriBead and Connector Ring)

### Características



**Filter**

Attach pluriStrainer® to a sterile 50 mL centrifuge tube. Then, add sample material onto the strainer and filter sample



**Reverse**

To obtain the larger fraction, take off the pluriStrainer®, turn it upside down onto another 50 mL tube and flush back the sample from the pluriStrainer®.



**Increase sample load**

With Funnel You can add up to 24 mL sample material on top.



**Stack**

Stacking of pluriStrainer® with different mesh sizes allows for straining various cell sizes at the same time.



**Flow control**

Allows to control the rate of flow by opening or closing the Luer-Lock, e.g. for physical dissociation of primary tissue (brain, spleen etc.).



**Low pressure**

If you add a syringe to the Connector Ring, it is possible to force low pressure to support the straining of rough sample material while pulling the piston.

**Additional Information**

Tiempo de entrega	1-2
El tamaño de la malla	100 µm
Embalaje	Individual single packed
Color	amarillo neón
Control de flujo	Yes, in combination with Connecting Ring
Mesh Fixation	Injected in housing for strong hold
Material de la tela	PET (Polietilentereftalato)

Material de la carcasa	LD-PE (Polietileno de baja densidad)
Mesh Type	woven
Manejo	Extended lip on strainer enables aseptic handling
Esterilidad	Estéril
Comparable con	Corning Cell Strainer, EASYstrainer™ Cell sieve, Falcon™ cell strainer, Reversible Strainer, MACS® SmartStrainer
Estabilidad	Buffer RLT®, TRIzol®, Isopropanol, organic solvent
Centrifugable	Sí
Disociación de tejidos	Sí
Desinfectable	Yes, with 70% ethanol
Área de interés	Inmunología, Biología celular, Biología vegetal, Microbiología, Oncología, Ecología
Condición de envío	Room Temperature
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,
Información legal	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center, Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning, Inc.; EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Falcon™ Cell Strainers is a registered trademark of Corning, Inc.; MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.