

## PLT Spin Medium (250 ml)

**Order No.:**60-00094-11



### Description

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

### Additional Information

|                      |  |
|----------------------|--|
| Size                 | 250 ml   |
| Delivery Time (days) | 1-2  |
| Storage Condition    | Room temperature, protect from light                         |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range             | 6.8 - 7.4  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.