

## universal S-pluriBead® anti-hamster

**Order No.:**31-GTaHS-10



### 描述

pluriBead universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

### Additional Information

大小	2 ml
品种	仓鼠
细胞类型	任何
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.