

## universal M-pluriBead® anti-rabbit

**Order No.:**31-GTaRB-20



### Descripción

pluriBead universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

### Additional Information

Tamaño	2 ml
Especies	Conejo
Tipos de células	Cualquier
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.