

CD193 S-pluriBead® anti-hu

Order No.:19-19300-10



Descripción

CD193 pluriBead® suspension for the positive cell isolation of CD193+ target cells from any human sample material. (2ml)

Additional Information

Tamaño	2 ml
Especies	Humano
Tipos de células	Granulocitos de basófilo, Granulocitos de eosinófilos
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.