

pluriSpin® Human B Cell Enrichment



Order No.	Size
19-06001-31	2 ml
19-06001-32	10 ml

Descripción

pluriSpin® Human B Cell Enrichment is designed to isolate pure, viable and untouched B cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, T cells, NK cells and granulocytes, while the negatively isolated human B cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified B cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Especies	Humano
Tipos de células	Células B
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>80%
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.