

Syringe Strainer 400 µm (12 pcs., non-sterile)

Order No.:43-71400-50



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

Mesh Size	400 µm
变化	12 pcs, non-sterile
大小	12 pcs.
流量控制	Yes
外壳材料	LD-PE (Low Density Polyethylen)
Mesh Type	woven
不育症	Non-Sterile
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.