

## CD326 S-pluriBead® anti-hu

**Order No.:**19-32600-10



### 概要

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any human sample material. (2ml)

### Additional Information

大きさ	2 ml
しゅるい	人間
細胞型	循環腫瘍細胞, CTC (循環腫瘍細胞, 上皮細胞)
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.