

## 10x PBS Solution (100 ml)

**Order No.:** 60-00010-10



### 상세 설명

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

## Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (10mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

Size	100 ml
Delivery Time (days)	1-2
Area of interest	Cell Biology
Sterility	Sterile
pH Range	pH 6.8
Shipping Condition	Room Temperature

**pluriSelect Life Science**

Deutscher Platz 5a, 04103 Leipzig, DE

Phone: +49 341 333 858 0 | Fax: +49 341 333 858 99

support@pluriselect.com

<https://www.pluriselect.com/>

Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.