

Funnel (50 pcs. , non-sterile)

Order No.:42-50000-50



Description

Funnel - Sterile device to support sample load onto pluriStrainer.

Additional Information

| | |
|----------------------|--|
| Size | 50 pcs. |
| Delivery Time (days) | 1-2 |
| Variation | 50 pcs, non-sterile |
| Shipping Condition | Room Temperature |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.