

pluriPlix



Order No.:50-01010-80

Description

Universal rotation mixer without own drive system for usage on top of a magnetic stirrer (not included in delivery). Particularly designed for the gentle mixing of cells. Allows for 2 mixing options: over-head and tilting-rotation (plus: adjustable rotation angle).

Additional Information

| | |
|----------------------|--|
| Size | 1 device |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.